



Handypette



*"We Produce Equipment
That Helps You
Produce Results"*



- Easy volume setting by simply turning the plunger. Does not snag gloves.
- Click Stop digital system for volume setting ensures no scales to align, so no unwanted volume alterations.
- High accuracy and precision guaranteed.
- Ergonomy, light weight and soft plunger movement provides comfort and minimizes hand fatigue.
- Two step plunger operation allows "reverse pipetting technique" useful for pipetting certain liquids.
- Tool kit provided as standard supply for in-house servicing and re-calibration.
- Long service assured from *Handypette* + Pipettes as the materials used in various parts offer durability and longlife.
- Each pipette has its individual Serial Number and a Calibration Report is enclosed with.
- Easy to calibrate and service.
- Specifications confirm to DIN 12650 & EN-ISO 8655 standards.

FULLY AUTOCLAVABLE

Accuracy at maximum



Handypette + Variable Volume Digital Micropipettes



Cat No.	Volume Range	Increment	Accuracy %	Precision (CV%)
700000	0.5 - 10 µl	0.1 µl	± 3.5 - 1.0	3.0 - 0.8
700020	2 - 20 µl	0.1 µl	± 3.0 - 1.0	1.5 - 0.4
700040	5 - 50 µl	1 µl	± 2.0 - 0.6	2.0 - 0.3
700060	10 - 100 µl	1 µl	± 3.0 - 0.8	1.0 - 0.2
700080	20 - 200 µl	1µl	± 1.8 - 0.6	0.7 - 0.2
700100	100 - 1000 µl	10 µl	± 0.9 - 0.8	0.5 - 0.2
700120	200 - 1000 µl	10 µl	± 1.0 - 0.8	0.5 - 0.2

Handypette + Fix Volume Digital Micropipettes

Cat No.	Volume Range	Accuracy %	Precision (CV%)
700180	1 µl	± 0.4	5
700200	2.5 µl	±2.4	2.5
700220	5 µl	± 1.5	1.4
700240	10 µl	± 0.9	0.8
700260	20 µl	± 0.6	0.5
700280	25 µl	± 0.6	0.5
700300	50 µl	± 0.6	0.4
700320	100 µl	± 0.5	0.3
700340	200 µl	± 0.8	0.3
700360	250 µl	± 0.8	0.3
700380	500 µl	± 0.8	0.3
700400	1000 µl	± 0.8	0.3

* Intermediate Fix volume pipettes may be custom-made to your needs.

FULLY AUTOCLAVABLE



Long term stability of calibration.
Re-calibration rarely but easily performed.



RV INSTRUMENTS Private Limited

visit us : www.rvinstruments.com

