



# Handypette<sup>+</sup> Multichannel

(8-Channel) Micropipette

*"We Produce Equipment  
That Helps You  
Produce Results"*



- High accuracy and precision guaranteed.
- Ergonomy, light weight and soft plunger movement provides comfort and minimizes hand fatigue.
- Easy volume setting by simply turning the plunger. Does not snag gloves.
- Click Stop digital system for volume setting ensures no scales to align, so no unwanted volume alterations.
- Two step plunger operation allows "reverse pipetting technique" useful for pipetting certain liquids.
- Long service assured from Handypette<sup>+</sup> Multichannel Pipettes as the materials used in various parts offer durability and longlife.
- Tool kit provided as standard supply for in-house re-calibration.
- Each pipette has its individual Serial Number and a Calibration Report is enclosed with.
- Easy to calibrate and service.
- Specifications confirm to DIN 12650 & EN-ISO 8655 standards.



## FULLY AUTOCLAVABLE



An ISO 9001 : 2000  
Certified Company  
Cer. No. RQ 912300


Accuracy at maximum

## TECHNICAL SPECIFICATION

Handypette<sup>+</sup> Digital Multichannel (8-Channel) Micropipette

Cat No.	Volume range	Increment	Accuracy%	Precision (CV%)
600100	0.5-10µl	0.1µl	± 3.5-1.0	3.0-0.8
600110	2-20µl	0.1µl	± 3.0-1.0	1.5-0.4
600120	5-50µl	1µl	± 2.0-0.6	2.0-0.3
600130	10-100µl	1µl	± 3.0-0.8	1.0-0.2
600140	20-200µl	1µl	± 1.8-0.6	0.7-0.2
600150	30-300µl	1µl	± 1.5-0.6	0.7-0.2
600160	50-300µl	1µl	± 1.5-0.6	0.8-0.2

### FULLY AUTOCLAVABLE



### Support Services for Better performance

RV Instruments (P) Ltd. provides complete maintenance services and all spare parts for the Digital Micropipette, has a wide network of distribution and sub-dealers which helps you getting better after sales services.



# RV INSTRUMENTS Private Limited

visit us : [www.rvinstruments.com](http://www.rvinstruments.com)



An ISO 9001 : 2000  
Certified Company  
Cer. No. RQ 912300

#### Warranty

All pipette manufactured by R V Instruments (P) Ltd. are warranted to be free from defects in material and workmanship for a period of one year, from the date of purchase. Routine cleaning and recalibration are not under the terms of warranty.