

# Mini Vertex Pipette Specification

## TECHNICAL SPECIFICATION

### Mini Vertex Fixed Volume Pipette

Cat No.	Volume Capacity	Accuracy %
900100	3 $\mu$ l	$\pm$ 2.4
900101	5 $\mu$ l	$\pm$ 1.5
900102	10 $\mu$ l	$\pm$ 0.9
900103	20 $\mu$ l	$\pm$ 0.6
900104	25 $\mu$ l	$\pm$ 0.6
900105	50 $\mu$ l	$\pm$ 0.6
900106	100 $\mu$ l	$\pm$ 0.5
900107	200 $\mu$ l	$\pm$ 0.8
900108	250 $\mu$ l	$\pm$ 0.8
900109	500 $\mu$ l	$\pm$ 0.8
9001010	1000 $\mu$ l	$\pm$ 0.8

### Mini Pipette Support



- The minipette is 85mm in length which is minimum size with maximum accuracy.
- Unique Tip-con design to accept tips up to 200ul.
- Two stage plunger movement ensuring last drop dispensing.
- Calibration through EN ISO 8655 standards.
- Full Accuracy, Quarter Size, Little Price